
PROVOX Adhesives

MD

CE

UK
CA

42°C
108°F
RT
2°C
36°F

1. Provox StabiliBase/StabiliBase OptiDerm


Figure 1.1


Figure 1.2


Figure 1.3


Figure 1.4


Figure 1.5


Figure 1.6


Figure 1.7


Figure 1.8

2. Provox FlexiDerm/XtraBase/OptiDerm


Figure 2.1


Figure 2.2


Figure 2.3


Figure 2.4

Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Provox Adhesives

Compatible Products; Kompatible Produkte; Compatibele producten; Produits compatibles; Prodotti compatibili; Productos compatibles; Produtos compatíveis; Kompatibla produkter; Kompatible produkter; Kompatible produkter; Yhteensopivat tuotteet; Samhæfar vörur; Ühilduvad tooted; Suderinami gaminiai; Kompatibilní produkty; Kompatibilis termékék; Kompatibilné produkty; Združljivi izdelki; Zgodne produkty; Produkse compatibile; Kompatibilni proizvodi; Συμβατά προϊόντα; Uyumlu Ürünler; Համատեղելի սարքերը; Совместимые изделия; Produk-Produk yang Kompatibel; Produk yang Sesuai; 互換性のある製品; 호환 가능한 제품; 相容產品; Sản phẩm tương thích; מוצרים תואמים; محمولات سازگار؛ تاقاوت مالا تاجت نام ل

	Provox XtraBase	Provox Stabilibase	Provox FlexiDerm	Provox OptiDerm
Provox XtraFlow HME	X	X	X	X
Provox XtraMoist HME	X	X	X	X
Provox ShowerAid	X	X	X	X
Provox Silicone Glue	X	X	X	X
Provox Adhesive Remover	X	X	X	X
Provox Skin Barrier	X	X	X	X
Provox Adhesive Strip	X	X	X	X
Provox Micron HME	X	X	X	X
Provox BasePlate Adaptor	X	X	X	X
Provox Cleaning Towel	X	X	X	X
Provox LaryTube with Ring	X	X	X	X

Contents

EN - ENGLISH.....	6	SK - SLOVENČINA.....	38
DE - DEUTSCH.....	8	SL - SLOVENŠČINA	40
NL - NEDERLANDS.....	10	PL - POLSKI.....	42
FR - FRANÇAIS.....	12	RO - ROMÂNĂ	44
IT - ITALIANO	14	HR - HRVATSKI.....	46
ES - ESPAÑOL.....	16	EL - ΕΛΛΗΝΙΚΑ	48
PT - PORTUGUÊS.....	18	TR - TÜRKÇE	50
SV - SVENSKA	20	HY - ՀԱՅԵՐԵՆ.....	52
DA - DANSK.....	22	RU - РУССКИЙ.....	54
NO - NORSK.....	24	ID - BAHASA INDONESIA ...	56
FI - SUOMI.....	26	MS - BAHASA MELAYU	58
IS - ÍSLENSKA	28	JA - 日本語.....	60
ET - EESTI.....	30	KO - 한국어.....	62
LT - LIETUVIŲ KALBA	32	ZHTW - 繁體中文.....	64
CS - ČESKY	34	TH - ไทย.....	66
HU - MAGYAR	36	VI - TIẾNG VIỆT.....	68
		HE - עברית.....	71
		AR - العربية.....	73
		FA - فارسی.....	75

1. Intended Use

Provox Adhesives are single-use devices intended for laryngectomized patients breathing through a tracheostoma. The devices are attached to the skin around the tracheostoma in order to provide attachment of components of the Provox HME System.

2. Description of the Device

Provox Adhesives are adhesives for everyday use; they are designed to ensure an airtight attachment for the Provox HME system components. The adhesives consist of an adapter and an adhesive part with peel-off liners.

Provox FlexiDerm is a very flexible material and has the strongest adhesive properties. It is a sticky, yet soft and flexible adhesive.

Provox OptiDerm is made of a hydrocolloid material that is suitable for sensitive skin. For example, dependent on clinical suitability and continued tolerance, it can be used postoperatively and during and after radiotherapy. The hydrocolloid material is commonly used in wound care and one of its main characteristics is that it absorbs wound fluids. When the material has absorbed fluids, it becomes white and, depending on the amount of secretions, it can become jelly-like, and it will no longer stick to the skin. Please note that the material may also absorb mucus coming from the stoma, and therefore it is still very important to always wipe and clear the mucus away from your adhesive and stoma area after you cough.

3. Warnings

- Do not use Provox FlexiDerm, Provox StabiliBase, or Provox XtraBase directly after laryngectomy.
- The use of Provox Adhesives needs to both be assessed on an individual basis and reviewed regularly through the course of radiotherapy treatment.
- Reuse will affect performance and may cause transfer of microorganisms leading to infections.
- Use only Provox Adhesives with compatible HMEs and accessories.

4. Precautions

- The adhesive may irritate the skin. Stop using the adhesive if skin irritation develops and consult your clinician.
- When cleaning the skin from residual glue, prevent particles/fluids from entering the tracheostoma.
- When using OptiDerm in the postoperative period or on sensitive skin, it should be removed very slowly and carefully.

5. How to Use

5.1 Prepare

1. Inspect that the adhesive is not damaged or has tear in the adhesive or around the base plate.
2. Always clean the skin with a Provox Cleaning Towel and/or soap and water.
3. Let the skin dry before applying the adhesive.

5.2 Apply

For guidance on how to apply, see Figures 1.1–1.8 and 2.1–2.4.

Apply skin protection products if needed, such as Provox Skin Barrier or adhesive supporting products such as Provox Silicone Glue (read the Instructions for Use accompanying the product).

5.3 Remove

1. Peel the adhesive off gently, using the finger lift tab. An adhesive remover, such as Provox Adhesive Remover (read the Instructions for Use accompanying the product), may be helpful for removing adhesives.
2. Always clean the skin with a Provox Cleaning Towel and/or soap and water after removal.
3. Dry the area carefully.

6. Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

7. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.
Atos Medical AB