
PROVOX Life

BasePlate Adaptor



MD


CE

UK CA


42 °C
108 °F
RT
2 °C
36 °F

ISO 15 mm


TrachPhone

50-60 °C
122-140 °F


1.

2 /250 ml (8,5 fl oz)
DISH SOAP
15 min
50-60 °C
122-140 °F


2.

3.

70% Ethanol 10 min
/(or)
70% Isopropylalcohol 10 min
/(or)
3% Hydrogenperoxide 60 min

4.

2 h

5.

Disclaimer

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Compatible Products; Kompatible Produkte; Compatibele producten; Produits compatibles; Prodotti compatibili; Productos compatibles; Produtos compatíveis; Kompatibla produkter; Kompatible produkter; Kompatible produkter; Yhteensopivat tuotteet; Samhæfar vörur; Ühilduvad tooted; Saderīgie izstrādājumi; Suderinami gaminiai; Kompatibilní produkty; Kompatibilis termékek; Združljivi izdelki; Kompatibilné produkty; Zgodne produkty; Produse compatibile; Kompatibilni proizvodi; Kompatibilni proizvodi; Συμβατά προϊόντα; Uyumlu Ürünler; Համատեղելի սարքեր; Совместимые изделия; Produk yang Sesuai; 互換性のある製品; 호환 가능한 제품; 相容產品; 兼容产品; מוצרים תואמים; المنتجات المتوافقة

Provox Life LaryTube
Provox Life LaryButton
Provox Life Standard Adhesive
Provox Life Sensitive Adhesive
Provox Life Night Adhesive
Provox Life Stability Adhesive
TrachPhone
Freevent DualCare
Freevent XtraCare
ISO 15 HMEs

Contents

EN - ENGLISH 7
DE - DEUTSCH 9
NL - NEDERLANDS 11
FR - FRANÇAIS 13
IT - ITALIANO 15
ES - ESPAÑOL 17
PT - PORTUGUÊS 19
SV - SVENSKA 21
DA - DANSK 23
NO - NORSK 25
FI - SUOMI 27
IS - ÍSLENSKA 29
ET - EESTI 31
LV - LATVIESU 33
LT - LIETUVIŲ KALBA 35
CS - ČESKY 37
HU - MAGYAR 39
SL - SLOVENŠČINA 41
SK - SLOVENČINA 43
PL - POLSKI 45
RO - ROMÂNĂ 47
HR - HRVATSKI 49
SR - SRPSKI 51
EL - ΕΛΛΗΝΙΚΑ 53
TR - TÜRKÇE 55
HY - ՀԱՅԵՐԵՆ 57
RU - РУССКИЙ 59
MS - BAHASA MELAYU 61
JA - 日本語 63
KO - 한국어 65
ZHTW - 繁體中文 67
ZHCN -简体中文 69
HE - עברית 72
AR - العربية 74

Intended use

Provox Life BasePlate Adaptor is an accessory that allows attaching medical devices, e.g. a HME, with an ISO 15 mm standard connector to a Provox Life attachment.

Contraindications

Shall not be used for mechanical ventilation.

Description of the device

Provox Life BasePlate Adaptor is an accessory that allows attaching and detaching medical devices, e.g. an HME, with an ISO 15 mm standard connector to a Provox Life attachment.

WARNINGS

- Unconscious patients must be constantly monitored. Attached devices may inadvertently slide off the adaptor, or the baseplate may detach from the tracheostoma.
- Make sure the adaptor is free from oil or any other lubricant before and during use. Lubrication may lead to detachment of devices from the adaptor, or of the adaptor from the baseplate.

PRECAUTIONS

- Read and follow the instructions of the medical device to be used with the adaptor.
- Never use force while attaching or detaching devices or the adaptor itself.
- Do NOT sterilize! The device is for single patient use only.
- Do not use if the device shows any signs of cracking or structural damage.

Cleaning and maintenance

Clean and disinfect the device before each use. Use drinking water to clean and rinse the device.

- Rinse the device in water (Fig. 1).
- Place the device in water-based mild soap solution for 15 minutes (Fig. 2).
- Clean inside of the device with Provox TubeBrush or Provox Swab. Rinse the device with water (Fig. 3).
- Disinfect the device at least once a day with one of the following methods: Ethanol 70% for 10 minutes or Isopropylalcohol 70% for 10 minutes or Hydrogenperoxide 3% for 60 minutes (Fig. 4).
- Let the device dry at least for 2 hours or until it looks dry (Fig. 5).

CAUTION: Do not use the device until it is completely dry. Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.
The device shall be cleaned and disinfected before each use to avoid risk of contamination. The device is for single patient use only, reuse between patients may cause cross contamination.
During hospitalization it is important to clean and disinfect the device immediately after use and again just before use due to an increased risk of infection and device contamination. Use sterile water for rinsing, rather than drinking water.

Device lifetime

Provox Life BasePlate Adaptor may be used for a maximum period of 3 months. If the device shows any signs of damage it shall be replaced earlier.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2026 Coloplast A/S. All rights reserved.
Atos Medical AB