Tracheofix

IFU

MD

Atos Medical AB
Kraftgatan 8,
SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00
info@atosmedical.com

Atos and the Atos logo are trademarks of Coloplast A/S.
© 2026 Coloplast A/S. All rights reserved.

Atos
atosmedical.com

LOT

Batch code; Chargencode; Partijnummer; Code de lot; Codice lotto; Código de lote; Código do lote; Batchkod; Batchkode; Eränumero; Šifra serije; Parti kodu

REF

Product reference number; Produktreferenznummer; Referentienummer product; Numéro de référence du produit; Numero riferimento prodotto; Número de referencia del producto;
Referência do produto; Referensnummer; Produktreferencenummer; Artikkelinumero;
Referentni broj proizvoda; Ürün referans numarası

Do not re-use; Nur zum Einmalgebrauch; Niet hergebruiken; Ne pas réutiliser; Non riutilizzare; No reutilizar; Não reutilizar; Återanvänd ej; Må ikke genanvendes; Ei uudelleenkäytettävä;
Nemojte ponovno upotrebljavati; Tekrar kullanmayın

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Valmistaja; Proizvođač; Üretici

DE

Country of manufacture with date of manufacture; Herstellungsland mit Herstellungsdatum; Land van productie met productiedatum; Pays de fabrication avec date de fabrication; Paese e data di produzione; País de fabricación con fecha de fabricación; País de fabrico com data de fabrico; Tillverkningsland med tillverkningsdatum; Fremstillingsland med fremstillingsdato; Valmistusmaa ja valmistuspäivämäärä; Država proizvodnje i datum proizvodnje; Üretim tarihi ile üretim ülkesi

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Anvendes inden; Viimeinen käyttöpäivä; Rok uporabe; Son kullanım tarihi

Keep away from sunlight and keep dry; Vor Sonnenlicht und Feuchtigkeit schützen; Niet blootstellen aan zonlicht, droog bewaren; Maintenir à l'abri de la lumière directe du soleil et de l'humidité; Tenere all'asciutto e al riparo dalla luce del sole; Mantener seco y alejado de la luz solar; Manter afastado da luz solar e da chuva; Håll borta från solljus och håll torr; Holdes bort fra sollys og holdes tør; Säilytettävä kuivassa paikassa suojassa auringonvalolta; Čuvajte na suhom mjestu i podalje od Sunčeve svjetlosti; Güneş ışığından uzakta ve kuru tutun

Storage temperature limit; Lagertemperaturgrenze; Maximale bewaartemperatuur; Limite de température de stockage; Temperature di stoccaggio; Límite de temperatura de almacenamiento; Limite da temperatura de armazenamento; Temperaturgränser för förvaring; Opbevaringstemperaturbegrænsning; Säilytyslämpötilan raja; Ograničenje temperature skladištenja; Saklama sıcaklığı sınırı

MAX
RT
MIN

Store at room temperature. Temporary deviations within the temperature range (max-min) are allowed; Bei Zimmertemperatur lagern. Temperaturschwankungen innerhalb des Temperaturbereichs (min.–max.) sind zulässig; Bij kamerte mperatuur bewaren. Tijdelijke afwijkingen van de temperatuur binnen de marge (min-max) zijn toegestaan; Stocker à température ambiante. Des écarts temporaires dans la plage de température ambiante (max-min) sont autorisés; Conservare a temperatura ambiente. Sono consentiti scostamenti temporanei entro i limiti di temperatura (max-min); Almacenar a temperatura ambiente. Se permiten desviaciones temporales dentro del rango de temperatura (máx.-mín.); Guardar à temperatura ambiente. São permitidos desvios temporários dentro do intervalo das temperaturas (máx-min); Förvara i rumstemperatur. Tillfälliga avvikelser inom temperaturområdet (max–min.) är tillåtna; Opbevar ved rumtemperatur. Midlertidige afvigelser i temperaturintervallet (maks-min) er tilladt; Säilytä huoneenlämmössä. Väliaikainen poikkeama lämpötila-alueella (suurin-pienin) on sallittu; Skladištite pri sobnoj temperaturi. Dopuštena su privremena odstupanja unutar raspona temperature (maks. – min.); Oda sıcaklığında saklayın. Sıcaklık aralığı (maksimum-minimum) içerisinde geçici sapmalara izin verilir

Caution, consult instructions for use; Vorsicht, Gebrauchsanweisung beachten; Let op: Raadpleeg de gebruiksaanwijzingen; Mise en garde, consulter le mode d'emploi; Attenzione, consultare le istruzioni per l'uso; Precaución, consultar instrucciones de uso; Cuidado, consultar as instruções de utilização; Varning, läs instruktion; Forsigtig, se brugsanvisning; Varoitus, katso käyttöohjeet; Pozor: pročitajte upute za upotrebu; Dikkat, kullanmadan önce talimatları inceleyin

MD

Medical Device; Medizinprodukt; Medisch hulpmiddel; Dispositif médical; Dispositivo medico; Producto sanitario; Dispositivo médico; Medicinteknisk produkt; Medicinsk udstyr; Lääkinnällinen laite; Medicinski proizvod; Tıbbi cihaz

IFU

Instructions for use; Gebrauchsanweisung; Gebruiksaanwijzing; Mode d'emploi; Istruzioni per l'uso; Instrucciones de uso; Instruções de utilização; Bruksanvisning; Brugsanvisning; Käyttöohjeet; Upute za upotrebu; Kullanım talimatları

Triman symbol and Infotri for France; Triman-Logo und Infotri für Frankreich; Triman- en Info-tri-symbool voor Frankrijk; Signalétique Triman et info-tri pour la France; Logo Triman e Infotri per la Francia; Símbolo Triman e Infotri para Francia; Símbolo Triman e Infotri para a França; Triman-symbol och Infotri för Frankrike; Triman-symbol og Infotri for Frankrig; Triman-symboli ja Infotri-ohjeet Ranskaa varten; Simbol Triman i oznaka Infotri za Francusku; Fransa için Triman sembolü ve Infotri

Recycling guidelines; Recycling-Richtlinien; Recycling richtlijnen; Consignes de recyclage; Linee guida per il riciclo; Directrices de reciclado; Orientações relativas à reciclagem; Riktlinjer för återvinning; Retningslinjer for genanvendelse; Kierrätysohjeet; Smernice za recikliranje; Geri dönüşüm yönergeleri;

CE

Indicates that the product is in compliance with European legislation for medical devices; Zeigt an, dass das Produkt mit den europäischen Rechtsvorschriften für Medizinprodukte konform ist; Geeft aan dat het product voldoet aan de Europese wetgeving voor medische hulpmiddelen; Indique que le produit est conforme à la législation européenne relative aux dispositifs médicaux; Indica che il prodotto è conforme alla legislazione europea sui dispositivi medici; Indica que el producto cumple la legislación europea sobre productos sanitarios; Indica que o produto está em conformidade com a legislação europeia relativa a dispositivos médicos; Anger att produkten överensstämmer med europeisk lagstiftning för medicintekniska produkter; Angiver, at produktet er i overensstemmelse med europæisk lovgivning for medicinsk udstyr; Osoittaa, että tuote on lääkinnällisiä laitteita koskevan EU-lainsäädännön mukainen; Označava da je proizvod u skladu s europskim zakonodavstvom za medicinske uređaje;

UK CA

UK Conformity Assessed (UKCA) marking; UK Conformity Assessed (UKCA)-Kennzeichnung; UKCA-markering (UK Conformity Assessed); Marquage UKCA (UK Conformity Assessed); Marcatura UK Conformity Assessed (UKCA); Marcado UKCA (UK Conformity Assessed);
Marcação (UKCA) de Conformidade Avaliada no Reino Unido; UKCA-märkning (UK Conformity Assessed); UKCA-mærkning (UK-overensstemmelsesvurdering); Yhdistyneen kuningaskunnan vaatimustenmukaisuuden arvioinnin (UKCA) merkintä; Oznaka ocjene sukladnosti sa zakonima UK-a (UKCA); Birleşik Krallık Uygunluk Değerlendirmesi (UKCA) işareti

UK Responsible Person
Coloplast Limited
Nene Hall, Peterborough Business Park
Peterborough
Cambridgeshire PE2 6FX

XXXXX, NN
YYYY-MM-DD

XXXX, NN = Reference number; Version number; Referenznummer; Versionsnummer; Referentienummer; Versienummer; Numéro de référence; Numéro de version; Numero di riferimento; numero di versione; Número de referencia; Número de versión; Número de referência; Número da versão; Referensnummer; Versionsnummer; Referencenummer; Versionsnummer; Viitenumero; Versionumero; Referentni broj; Broj verzije; Referans numarası; Sürüm numarası

YYYY-MM-DD = Date of issue; Ausstellungsdatum; datum van uitgifte; Date d'émission; Data di emissione; Fecha de expedición; Data de emissão; Utfärdandedatum; Udstedelsesdato; julkaisupäivä; Datum izdaje; Verilme tarihi;

EN - ENGLISH

Intended Use

Provox Tracheofix is a single use foam protector intended to absorb secretions and to provide protection and aesthetic coverage of the tracheostoma.

Instructions for Use

Single use only: do not reuse, this may increase the risk of possible transmission of infection.
Sensitive skin can increase the risk of skin irritation/breached skin.
Make sure the stoma area is clean and dry before application.
Remove the paper strips from the device and apply the device to the skin. See Fig. 1 and Fig. 2.
Change the device at least every 24 hours, or more often if needed. Be careful when removing the device to not cause skin damage.

Information

The products are intended to be stored in the box until removed for use as the box protects the products, and the label on the box contains important information like batch code, used-by date, and a unique barcode.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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